

colorimeter experiment of Cu^{2+} ☆ ☰ ☰

File Edit View Insert Format Tools Add-ons Help Last edit was on July 5, 2021

100% Normal text Arial 11 B I U A

Google Account
Deepak Nagare
deepaknagare03@gmail.com

Observation Table :-
Determination of λ_{max} of given solution of Cu^{2+} ions and verify Beer's law.

Filter No.	Wavelength λ_{nm}	Absorbance (A)
1	400	0.12
2	420	0.17
3	440	0.19
4	460	0.24
5	480	0.19
6	500	0.17
7	520	0.11

Determination of Absorbance at different conc. for λ_{max}

Sr. No.	Concentration (M)	Absorbance or O.D. (A)	% Transmittance (T)
1	0.001	0.11	75.1
2	0.002	0.17	62.1
3	0.003	0.20	47.1
4	0.004	0.26	40.1
5	0.005	0.31	39.1
6	0.006	0.39	34.1
7	0.007	0.47	33.1
8	Unknown	0.53	43.1

Preparation of solutions :-

Stock solution of 0.01 M CuSO_4 is supplied to you. Using this solution prepare solutions of different concentration in volumetric flask as given in table.

Sl. No.	Concentration (M)	ml of 0.01 M CuSO_4	ml of Dist. Water	Total Volume (ml)
1	0.001	10	90	100
2	0.002	20	80	100
3	0.003	30	70	100
4	0.004	40	60	100
5	0.005	50	50	100
6	0.006	60	40	100
7	0.007	70	30	100
8	Unknown	Unknown	20	100

https://accounts.google.com/SignOutOptions?hl=en&continue=https://docs.google.com/document/d/1XQzHob8kx2pM_2AJSld-eL3bK8YgTos7CveeHVM/edit#3Foid%3D107078427122713791543%26usp%3Ddocs_home%26th%3Dtrue&service...